

JOBST® Elvarex®



**Custom-made
garments for the
management of
Lymphedema**

15-90 mmHg*

#1 physician recommended

JOBST®

Comfort, Health, and Style!



Elvarex® features and benefits:

JOBST® Elvarex® helps maintain the reduced volume and shape of the limb. These garments help to inhibit further fluid accumulation and help to assure therapeutic results achieved by compression. Elvarex compression garments are indicated when decongestion is achieved so the garment helps prevent an increase in the extremity dimensions.

- Custom manufacturing to ensure a contoured fit to the shape of the limb
- Controlled gradient compression for the effective long-term management of lymphedema
- Unique knitted construction that allows greater air permeability for superior comfort and easier garment care
- Excellent durability for extended wear
- Upper extremity garments in class 1, 2 and 2 Forte
- Lower extremity garments in class 1, 2, 3, 3 Forte, 4 and 4 Super
- Flat knitting process with seam provides superior fit, especially for elbow region in sleeves
- Diagonal cut on top of sleeve
- Slip form cut on lower extremity garments
- Slant open toe
- Available in a variety of styles and options

Available colors:



Actual garment color may vary from printed swatch sample.

*Upper extremity ONLY

Indications:









- Management of lymphedema and other edema

Contraindications:









- Ischemia (e.g. advanced arterial disease)
- Uncontrolled congestive heart failure
- Untreated septic phlebitis
- Phlegmasia cerulea dolens
- Dermatitis in the acute phase
- Arterial insufficiency
- Wet dermatosis
- Cutaneous infections

Elvarex styles and options:








Upper extremity styles:

							
Sleeve, wrist to axilla	Sleeve with shoulder cap and strap	Sleeve with shoulder cap and bra loop, with Velcro®	Gauntlet with thumb and without fingers	Glove with thumb and fingers	Sleeve with gauntlet/glove attached	Sleeve with shoulder cap and strap, with gauntlet/glove attached	Sleeve with shoulder cap and bra loop with Velcro, with gauntlet/glove attached

Lower extremity styles:

							
Knee (AD) with slant open toe	Thigh (AG) with slant open toe and slip form	Chap style (AGT) slant open toe with slip form	One-Leg Pantyhose (AT) consists of a waist-high pantyhose with 1 1/2 leg	Two-Leg Pantyhose (AT)	Two-Leg Pantyhose (AT) with fly	Capri	Bermuda

Options

						
Two-Leg Pantyhose with adjustable waist band	Foot Cap	Slanted, Closed Toe	Straight, Open Toe	Straight, Closed Toe	Silicone Band (available on upper extremity sleeves, knee and thigh)	Zippers (available on select garments)

Please refer to the product label and package insert for full instructions on the safe use of this product.

CAUTION: This product contains natural rubber latex which may cause allergic reactions.

CAUTION: Federal U.S.A. Law restricts the sale of this device to or on the order of a physician. Velcro® is a registered trademark of Velcro U.S.A., Inc.

* = Design Pressure

www.jobst.com



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.jobst.com
To order toll-free: JOBST (+1) 800 537 1063

50838 R7

© 2013 BSN medical Inc.

 physician recommended



Comfort, Health and Style!

REV 01/13