

JOBST Confidence® Lower Extremity Order Form

TO ORDER:
hms-elvarex-orders@essity.com
<https://eshop.jobst-usa.com>
Tel: (+1) 800-537-1063
Fax: (+1) 800-835-4325

Patient Name / ID Code or File # _____ Date _____
 Address _____ Gender M F
 City/State/Zip Code _____ DOB _____

Fitter Name _____	Fitter # _____ <small>(Not Required)</small>	Fitter Phone _____
Fitter Facility _____	Email _____	
Ship To Acct # _____	Acct Name _____	
Address _____	City _____	State _____ Zip Code _____
Email _____	Phone _____	Fax _____
Bill To Acct # _____	Acct Name _____	
Address _____	City _____	State _____ Zip Code _____
Email _____	Phone _____	Fax _____

PO# _____

Original Order Reorder w Changes
 Exact Reorder Schema # _____

Confirmation Fax # _____

Email _____

By choosing communication via email (above), I acknowledge that Personal Health Information associated with this purchase may be transmitted from Essity in a non-encrypted manner.

Color	Styles	Quantity/Class	CCL1 <small>(18-21mmHg†)</small>	CCL2 <small>(23-32mmHg†)</small>	CCL3 <small>(34-46mmHg†)</small>
<input type="checkbox"/> Beige <input type="checkbox"/> Anthracite Heather <input type="checkbox"/> Black <input type="checkbox"/> Jeans Heather <input type="checkbox"/> Caramel <input type="checkbox"/> Red Heather	<input type="checkbox"/> AD Knee <input type="checkbox"/> AB1 Sock <input type="checkbox"/> CT Capri <input type="checkbox"/> AG Thigh <input type="checkbox"/> BT Capri <input type="checkbox"/> ET Bermuda <input type="checkbox"/> AT Panty <input type="checkbox"/> B1T Capri <input type="checkbox"/> AG-HT 1 Leg Panty	Left _____ Right _____			

Straight Open Toe Length††
 Lateral IA _____ cm

Straight Closed Toe Length††
 Total Foot IZ _____ cm

Slant Open Toe Length
 Medial IA _____ cm
 Lateral IA _____ cm

Slant Closed Toe Length
 Medial IA _____ cm
 Lateral IA _____ cm
 Total Foot IZ _____ cm


AD Band Options
 Without Silicone
 SoftFit Band AD
NOTE: this is a 5cm band

AG Band Options
 5cm Dotted Band With Lateral Rise (Standard)

Measuring Guidelines
(Only applicable for Confidence) See box below for applicable tension at each landmark.

0 no tension
 + light tension
 ++ heavy tension

Options	Circum. (c)		Length (l)	
	Left	Right	Left	Right
<input type="checkbox"/> Lateral Rise =10% of circumference at D/G and is not adjustable (ex: if cD/ cG is 35cm then lateral rise is 3.5cm)*			/TT	/T
<input type="checkbox"/> Men's style <input type="checkbox"/> with fly <input type="checkbox"/> without fly			/K3	/H
<input type="checkbox"/> Floral Waistband Women 5cm			/K	
<input type="checkbox"/> Elastic Waistband Women 5cm			/G	
<input type="checkbox"/> Elastic Waistband Men 4cm			/F	
<input type="checkbox"/> Decorative Line (Front of garment)			/E1	
<input type="checkbox"/> Patient Initials Max 2 letters (A-Z)			/E	
<input type="checkbox"/> Ankle Comfort Zone			/D	
<input type="checkbox"/> Knee Comfort Zone			/C	
<input type="checkbox"/> Hallux Valgus (slant toe option only)			/B1	
			/B	
			/Y	
			/B1++	
			/B	
			/Y	
			/A++	



* only available in AD and AG
 † design pressure
 Δ cE1 for Bermuda only, measure 4cm above kneecap
 ** silicone band & straight ending
 †† n/a Hallux Valgus
 *** if silicone band
 ΔΔ if patient is standing

essity BSN Medical Inc., an Essity company
 5825 Carnegie Blvd. Charlotte, NC 28209-4633

Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application

Size	Shoe size	BNR	UOM / Box	Order Qty
X-Small	≤ 2	7965803	1	
Small	2.5 - 5.5	7965804	1	
Medium	6 - 8	7965802	1	
Large	8.5 - 11	7965902	1	
X-Large	≥ 11.5	7966001	1	



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	≤ 7.5	7957100	1	
Large	8 - 12	7957101	1	
X-Large*	≥ 13	7957102	1	



Arion Magnide® 2in1

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	< 7.5	7957110	1	
Large	≥ 8 - 10	7957111	1	
X-Large*	≥ 10.5	7957112	1	



* XL also suitable for extra wide feet with cA up to 28 cm

Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product

Size	Shoe size	BNR	UOM / Box	Order Qty
Small	< 3 - 6	7063301	1	
Medium	6.5 - 8	7063302	1	
Large	8.5 - 10.5	7063303	1	
X-Large	≥ 11	7063304	1	



Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children

Size	Age	Length	BNR	UOM / Box	Order Qty
X-Small	0-3 yrs	11.6" (29.5cm)	7957108	1	
Small	3-6 yrs	15.5" (39.5cm)	7957109	1	




Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®

Size	Circumference of lower leg	BNR	UOM / Box	Order Qty
Medium	≤ 15.7" (40cm)	7957106	1	
Large	≤ 23.6" (60cm)	7957107	1	
Extra large	≤ 39.5" (100cm)	7957107	1	



Dycem® Non-Slip Material*

- Improves grip
- Provides stability
- Helps with easier donning of a JOBST® garment

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
One size	7763600	1	



Dycem® is a registered trademark of Dycem Ltd.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
X-Small	7740027	1 pair	
Small	7740028	1 pair	
Medium	7740029	1 pair	
Large	7740030	1 pair	
X-Large	7740031	1 pair	



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box	Order Qty
Standard	110913	1 donner	
Large	110910	1 donner	

