



Confidence® Order Form Armsleeve

TO ORDER:

<https://eshop.jobst-usa.com/>

Tel: (+1) 800-537-1063 Fax: (+1) 800-835-4325

Patient Name / ID Code or File # _____ DOB _____ Date _____

Address _____ Gender M F

City/State/Zip _____

Diagnosis _____

Doctor/Address _____

City/State/Zip _____

PO#	
Original Order <input type="checkbox"/>	Reorder w Changes <input type="checkbox"/>
Exact Reorder <input type="checkbox"/>	Schema # _____

Fitter Name _____ Fitter # _____ Fitter Phone _____

Fitter Facility _____ Email _____

Ship To Acct # _____ Acct Name _____

Address _____ City _____ State _____ Zip _____

Email _____ Phone _____ Fax _____

Bill To Acct # _____ Acct Name _____

Address _____ City _____ State _____ Zip _____

Email _____ Phone _____ Fax _____

Confirmation Fax # _____

Email _____

By choosing communication via email (above), I acknowledge that Personal Health Information associated with this purchase may be transmitted from BSN in a non-encrypted manner.

Last 4 digits of credit card on file OR Exp. _____

New card - call to provide credit card # Billing Zip _____

Name on CC _____

Color

Beige Caramel
 Black Jeans Heather
 Red Heather Anthracite Heather

Styles

CG1
 AG1

Quantity/Class

	CCL1 (15-21mmHg*)	CCL2 (23-32 mmHg*)
Left		
Right		

Lateral Rise is 10% of circumference at G and is not adjustable (ex: if cG is 35cm then lateral rise is 3.5cm)

Elbow Options

Elbow Comfort Zone

Elbow Bend Options

Elbow 25 Degree (standard)
 Elbow 45 Degree

Decorative Options


Decorative Line (Front of garment)
 Patient Initials
 Max 2 letters (A-Z) _____

Silicone Band

No Silicone
 SoftFit 2.5cm
 2.5cm Top
 2.5cm Inside 1/2

Circumference Measurements (cms)


Length Measurements (cms)



Measuring Guidelines

(Only applicable for Confidence)
 See Arm Diagram for applicable tension at each landmark.
 Apply lay-on measurement/no tension to every circumference on the hand, fingers and thumb. See JOBST Measuring Guidelines for complete measuring instructions.

- 0 no tension
- + light tension
- cG = 0 no tension with silicone band
- cG = + light tension without band
- lC1 = 5 to 7cm above cC
- (lCG must be taken with the arm bent)



* Design Pressure
 For additional product order forms, please go to:
<http://www.jobstcompressioninstitute.com/resources/orders>



BSN Medical Inc., an Essity company
 5825 Carnegie Blvd. Charlotte, NC 28209-4633