



# Classic Arm Sleeve AG1 Ready-to-Wear

## TO ORDER:

Email:

ca.customerservice@essity.com

Tel: 1-877-978-5526

Fax: 1-877-978-9703

Patient Name: \_\_\_\_\_

Previous Patient?  Yes Gender:  F  M

Height\*: \_\_\_\_\_ Weight\*: \_\_\_\_\_ Birthdate: \_\_\_\_\_

All JoViPak RTW garments are pre-made for secondary lymphedema

\*Height and weight are required.

### PAYMENT INFORMATION

Account # <input type="checkbox"/> Bill to Account	Date
<input type="checkbox"/> Charge Credit Card <input type="text"/> <input type="text"/> Card Exp. Date	PO #
Card #	Fax Confirmation #
Name on Card	Email Confirmation

### BILLING ADDRESS

Business Name	Address	Attention	City	Province	Phone	Postal Code
---------------	---------	-----------	------	----------	-------	-------------


### SHIPPING ADDRESS

Same as Billing Address

Business Name	Address	Attention	City	Province	Phone	Postal Code
---------------	---------	-----------	------	----------	-------	-------------

### ORDER SPECIFICATIONS

Quote Only  Quote & Proceed



<b>Polartec® Power Dry® Colours &amp; Channeling</b>
<input type="checkbox"/> Royal Blue - Standard Channeling

<b>JoViJacket</b>
<input type="checkbox"/> Black Nylon & Spandex Powernet

(JoViJackets are required to be worn with your JoVi foam garment to ensure maximum fit and effectiveness)

<b>Side</b>
<input type="checkbox"/> Right Arm
<input type="checkbox"/> Left Arm

<b>Size</b>
<input type="checkbox"/> Small
<input type="checkbox"/> Medium
<input type="checkbox"/> Large
<input type="checkbox"/> XLarge

<b>Length</b>
<input type="checkbox"/> Regular
<input type="checkbox"/> Long

<b>Additional Items</b>
<input type="checkbox"/> Dycem® - donning aid
<input type="checkbox"/> Arion Easy Slide - donning aid

Dycem® is a registered trademark of Dycem Ltd.

Comments/Quantity:

Fitter/Therapist Name: \_\_\_\_\_ Phone: \_\_\_\_\_ Email: \_\_\_\_\_



JOBST®, an Essity brand

jobstcanada.com

Essity  
1275 North Service Road West, Suite 800  
Oakville, ON Canada L6M 3G4  
Tel. 1-877-978-5526 Fax 1-877-978-9703

63962E R1 136459 ©Essity Hygiene and Health AB, 2022 122

# Classic Arm Sleeve AG1

## Ready-to-Wear

Patient Name: \_\_\_\_\_

All measurements are in centimeters (cm).

Classic Arm Sleeve AG1 Size Chart			
	Least Wrist Circumference	Least Elbow Circumference	Axilla Circumference
	"C" Landmark	"E" Landmark	"G" Landmark
Small	15.5–16.5	23–26	27.5–32
Medium	16.5–18	26–29	32–36
Large	18–19	29–33	36–41
XLarge	19–20.5	33–37.5	41–45

Classic Arm Sleeve AG1 Length Chart	
Wrist to Axilla (C to G) (measured with arm hanging down at side)	
Regular	43
Long	45.5



### Arion Easy-Slide Arm

- The user-friendly application aid makes putting on compression arm sleeves quick and easy
- A straightforward donning method in combination with the application aid



Size	Circumference of widest part of the arm	BNR	UOM / Box	Order Qty.
Medium	37–38.5cm (14.5"–15.1")	7966102	1	
Large	39–41cm (15.3"–16.1")	7510001	1	