



Confidence® Lower Extremity Order Form

Patient Name / ID Code or File # _____ DOB _____ Date _____
 Address _____ Gender M F
 City/State/Zip Code _____
 Diagnosis _____ Lymphedema Lipedema

TO ORDER:
 hms-elvarex-orders@essity.com
 https://eshop.jobst-usa.com
 Tel: (+1) 800-537-1063
 Fax: (+1) 800-835-4325

Fitter Name _____ Fitter # _____ Fitter Phone _____
 (Not Required)
 Fitter Facility _____ Email _____
 Ship To Acct # _____ Acct Name _____
 Address _____ City _____ State _____ Zip Code _____
 Email _____ Phone _____ Fax _____
 Bill To Acct # _____ Acct Name _____
 Address _____ City _____ State _____ Zip Code _____
 Email _____ Phone _____ Fax _____

PO# _____
 Original Order Reorder w Changes
 Exact Reorder Schema # _____
 Confirmation Fax # _____
 Email _____
 By choosing communication via email (above), I acknowledge that Personal Health Information associated with this purchase may be transmitted from Essity in a non-encrypted manner.

Color	Styles	Quantity/Class	CCL1 (18-21mmHg†)	CCL2 (23-32mmHg†)	CCL3 (34-46mmHg†)
<input type="checkbox"/> Beige <input type="checkbox"/> Anthracite Heather <input type="checkbox"/> Black <input type="checkbox"/> Jeans Heather <input type="checkbox"/> Caramel <input type="checkbox"/> Red Heather	<input type="checkbox"/> AD Knee <input type="checkbox"/> AB1 Sock <input type="checkbox"/> AG High <input type="checkbox"/> BT Capri <input type="checkbox"/> AT Panty <input type="checkbox"/> B1T Capri	<input type="checkbox"/> ET Bermuda <input type="checkbox"/> AG-HT One Leg Panty			
		Left			
		Right			

AD Band Options
 Without Silicone
 SoftFit Band AD
 NOTE: this is a 5cm band

AG Band Options
 5cm Dotted Band With Lateral Rise (Standard)


Measuring Guidelines
 (Only applicable for Confidence) See box below for applicable tension at each landmark.

0 no tension
 + light tension
 ++ heavy tension

Options

- Lateral Rise =10% of circumference at D/G and is not adjustable (ex: if cD/ cG is 35cm then lateral rise is 3.5cm)*
- Men's style with fly without fly
- Floral Waistband Women 5cm
- Elastic Waistband Women 5cm
- Elastic Waistband Men 4cm
- Decorative Line (Front of garment)
- Patient Initials Max 2 letters (A-Z)
- _____
- Ankle Comfort Zone
- Knee Comfort Zone
- Hallux Valgus (slant toe option only)

	Circum. (c)	Length (l)	Length (l)	
cT ⁰		/TT	/T	
cH ⁰		/K3	/H	
	Circumference (c)		Length (l): Taken from each landmark to floor	
	Left	Right	Left	Right
			/K	
cG ^{++/+**}			/G	
cF ⁺⁺			/F	
cE1 ^Δ			/E1	
cE ⁺			/E	
cD ^{+/0...}			/D	
cC ⁺⁺			/C	
cB1 ⁺⁺			/B1	
cB ⁺			/B	
cY ⁰				
cA ^{+ΔΔ}				




* only available in AD and AG
 † design pressure
 Δ cE1 for Bermuda only, measure 4cm above kneecap
 ** silicone band & straight ending
 †† n/a Hallux Valgus
 *** if silicone band
 ΔΔ if patient is in lying position

essity
 BSN Medical Inc., an Essity company
 5825 Carnegie Blvd. Charlotte, NC 28209-4633

Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application


Size	Shoe size	BNR	UOM / Box	Order Qty
X-Small	≤ 2	7965803	1	
Small	2.5 - 5.5	7965804	1	
Medium	6 - 8	7965802	1	
Large	8.5 - 11	7965902	1	
X-Large	≥ 11.5	7966001	1	



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes



Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	≤ 7.5	7957100	1	
Large	8 - 12	7957101	1	
X-Large*	≥ 13	7957102	1	



Arion Magnide® 2in1

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	< 7.5	7957110	1	
Large	≥ 8 - 10	7957111	1	
X-Large*	≥ 10.5	7957112	1	





* XL also suitable for extra wide feet with cA up to 28 cm

Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product


Size	Shoe size	BNR	UOM / Box	Order Qty
Small	< 3 - 6	7063301	1	
Medium	6.5 - 8	7063302	1	
Large	8.5 - 10.5	7063303	1	
X-Large	≥ 11	7063304	1	



Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children

Size	Age	Length	BNR	UOM / Box	Order Qty
X-Small	0-3 yrs	11.6" (29.5cm)	7957108	1	
Small	3-6 yrs	15.5" (39.5cm)	7957109	1	




Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®

Size	Circumference of lower leg	BNR	UOM / Box	Order Qty
Medium	≤ 15.7" (40cm)	7957106	1	
Large	≤ 23.6" (60cm)	7957107	1	
Extra large	≤ 39.5" (100cm)	7957107	1	



Dycem® Non-Slip Material*

- Improves grip
- Provides stability
- Helps with easier donning of a JOBST® garment

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
One size	7763600	1	



Dycem® is a registered trademark of Dycem Ltd.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
Small	131202	1 pair	
Medium	131203	1 pair	



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box	Order Qty
Standard	110913	1 donner	
Large	110910	1 donner	

