



JOBST® JoViPak® Garments

These handmade compression garments are indicated for lymphedema patients (primary or secondary lymphedema) and other forms of edema who want a more comfortable and easier alternative to bandaging with complementary compression to daily recommended lymphedema therapy, during rest and night. The continuous flow channeling and blended density chip foam construction are designed to address all levels of lymphedema, including fibrotic tissue.

JOBST® JoViPak® in combination with the JOBST® JoViJacket is designed for lymphedema compression therapy during the night. The 4-blend foam mix and stitched channels are designed to create multiple pressure points and stimulate lymphatic drainage. These garments for edema management come in various styles and, with the addition of specialized pads, are designed to provide comfort and control from head to toe.

The **JOBST® JoViPak® Lower Leg Garments** and **JOBST® JoViJackets** are for patients who have swelling in their lower leg(s) and feet from lymphedema and chronic venous insufficiency. The sole of the foot extends to the tip of the longest toe, with the dorsum cut back at an angle to expose the toes and ends at the distal knee. The foot length is proportional to the leg length and girth.

The Ready-Made Classic Lower Leg Garment (AD) can also be used with a Hip Hugger or Boxer Capri. This combination offers great flexibility, particularly when there is a significant disparity between upper and lower leg sizes.

Available as Ready-Made (RTW) and Made-Measure (MTM)



MTM AG1

RTW AD

The **JoViPak® Hip Hugger** is intended for patients with thigh swelling extending up into the hip or with heavy thigh lobules that need lifting. It extends from just below the knee to the waist and is held in position by the unpadded Nylon/Spandex® wrap-around abdominal stabilizer, which wraps around the hip and is secured in front with Velcro®.

Because the Hip Hugger has an attached Velcro® wrap-around band, the garments are only available in Organic Cotton/Lycra® with a Nylon/Spandex® band. This garment can be modified to accommodate both legs. A good choice for patients desiring an open crotch, with no genital swelling, for toileting issues.

Hip Huggers are usually worn in combination with JoViPak® lower leg garments and full-length Nylon/Spandex® JoViJackets. The unpadded cuff on the lower end of the Hip Hugger fits inside the lower leg garment.

Available as Made-To-Measure (MTM)



The **JoViPak® Boxer** extends from mid-thigh to high waist (transverse watershed) and is available for both males and females. Ideal for the patient with hip, gluteal, and genital swelling. Custom Boxers end at mid-thigh or approximately 6" below groin, unless otherwise specified. All boxers come with zippered opening. If the Boxer will be worn with a leg garment (AF1 w/cuff), a full length JoViJacket can be ordered and will help hold the two garments together as one unit. Refer to measuring form and note that longer lengths require additional measurements.

The **JoViPak® Boxer Capri** can be ordered if swelling extends from below the knee into the hip up to the waist. It has a cuff at the "D" measurement so it can be worn with a Classic Lower Leg garment. This style may be more appropriate for patients with lobules or when the thigh area is more affected.

Both styles must be worn with a JoViJacket for optimal therapeutic benefit.

Available as Made-To-Measure (MTM)



Boxer

Boxer Capri

JOBST® JoViPak® Garments

The **JoViPak® Vest** is a versatile garment that addresses chest wall edema and fibrosis common to Lumpectomy and Mastectomy patients. The Vest has adjustable shoulder straps and is made of Cotton/Lycra® for comfort, allowing for chest expansion during normal respiration. The Velcro® midline closure separates at the waistline to eliminate buckling. A Nylon/Spandex® JoViJacket is recommended for optimal therapeutic benefit. It also has adjustable shoulder straps and wide adjustable Velcro® opening, so garment can be tightened down to patient's comfort level.

Available as Made-To-Measure (MTM)



The **JoViPak® Shoulder-Torso** garment is an innovative solution for pain and swelling following breast cancer surgery when there is metastasis into the brachial plexus. The design evolved to offer a separate zippered handpiece, a lower arm zipper for ease of donning, and an attached sling to support the dependent limb.

Available as Made-To-Measure (MTM)



The **JoViPak® Busti** is ideal for the lumpectomy patient who has swelling in the affected breast. A modification of the JoViPak® Vest, the Busti has no shoulder straps and extends down to the lowest rib, rather than to the waist.

Not recommended for over DD/E sized women. Mastectomy patients are better served by a Vest, Double Mastectomy pad, and/or Bellisse®. An outer Nylon/Spandex® JoViJacket will enhance the effectiveness of the Busti and secure it in position and is highly recommended when ordering this garment.

Available as Made-To-Measure (MTM)



The **JOBST® JoViPak® Arm** and **JOBST® JoViJacket** are designed to eliminate the need for nighttime bandaging of post-operative arm swelling for mastectomy and lumpectomy patients. The arm sleeve is contoured to cover the DIP's (distal portion of the fingers) and thumb to the base of the nail bed and features an anatomically positioned thumb stub and an antecubital dart for added comfort. Continuous Flow™ channeling directs lymph flow laterally, away from the axilla to functional alternate pathways.

Available as Ready-Made (RTW) and Made-To-Measure (MTM)



Micro-massage: Multiple pressure points designed to stimulate lymph flow and soften fibrotic tissue.

Massage effect via multiple densities of chipped foam