




# Flat-Knit Decision Table:

JOBST® Confidence vs JOBST® Elvarex Soft vs JOBST® Elvarex

Body Shape			
Stage	Mild (I)	Moderate (II)	Severe (III)
Patient Description	First time user, newly diagnosed, preventative treatment, Lipedema/ Lymphedema	First time user, low tolerance to stiff fabric, desire to “feel” the compression, Lipedema/ Lymphedema	Persistent/aggressive edema, Lymphedema (not best for Lipedema)
Skin Folds	None	Minimal	Significant
Tissue	Proportionate, symmetrical, athletic, controlled/maintained edema	Adipose present, beginning development of fibrosis	Significant adipose and/or hard, fibrotic
User Need	Light & Soft material	Soft & Flexible material	Firm & Strong material
JOBST Product Recommendation	<div>JOBST Confidence</div> <div>JOBST Elvarex Soft</div> <div>JOBST Elvarex</div>		