

Elvarex® Soft Order Form

Lower Extremity

TO ORDER:
<https://eshop.jobst-usa.com/>
Tel: (+1) 800-537-1063
Fax: (+1) 800-835-4325


Patient Name / Essity File # _____ DOB _____
 Address _____ Gender M F
 City/State/Zip _____
 Diagnosis _____ Date _____

Color	Quantity/Class	CCL1 18-21 mmHg*	CCL2 23-32 mmHg*	CCL3 34-46 mmHg*
<input type="checkbox"/> Beige <input type="checkbox"/> Gray	Left			
<input type="checkbox"/> Black <input type="checkbox"/> Cocoa	Right			
<input type="checkbox"/> Cranberry <input type="checkbox"/> Cherry	Body Bandage <small>CCL must be same as legs</small>			
<input type="checkbox"/> Navy				

Styles <input type="checkbox"/> AD Knee <input type="checkbox"/> AG-T Chap: <input type="checkbox"/> pc. <input type="checkbox"/> pr. <input type="checkbox"/> AG Thigh <input type="checkbox"/> AT Pantyhose <small>AT Pantyhose must be all one compression class. All leg lengths must be equal.</small>		<input type="checkbox"/> Straight Open Toe Length Lateral _____ cm <input type="checkbox"/> Straight Closed Toe Length Total Foot /Z _____ cm	<input type="checkbox"/> Slant Open Toe Length Medial _____ cm Lateral _____ cm	<input type="checkbox"/> Slant Closed Toe Length Medial _____ cm Lateral _____ cm Total Foot /Z _____ cm
---	--	--	--	--

Circum. (c)	Length (l)	Length (l)	
cT ⁰	/GT	/T	
cH ⁺⁺	/KT	/H	
Circumference (c)		Length (l): Taken from each landmark to floor	
Left	Right	Left	Right
		/K	
cG ⁺⁺		/G	
cF ⁺⁺		/F	
cE ⁺		/E	
cD ⁺⁺⁺		/D	
cC ⁺⁺		/C	
cB1 ⁺⁺		/B1	
cB ⁺		/B	
cY ⁰		/A (medial)	
cA ⁺		/A (lateral)	

Variations	Tensions
<input type="checkbox"/> B1G-T <input type="checkbox"/> BG-T <input type="checkbox"/> FT Biker Short	0 no tension + light tension ++ heavy tension



Special Options	
<input type="checkbox"/> Oblique: AD standard 4cm	<input type="checkbox"/> Adj. waistband (AT panty only)
<input type="checkbox"/> Other: ___cm (2cm-6cm) (AG fixed 7cm, no modifications)	<input type="checkbox"/> Open pubis (AT panty only)
	<input type="checkbox"/> T-Heel

Silicone Band	On Top
2.5cm (A-D Only)	
5cm	
AG-T Not available with Silicone band.	
AT Pantyhose must be all one compression class. All leg lengths must be equal.	
<input type="checkbox"/> SoftFit band (A-D Only)	
Pocket	Lining (Pocket all sides closed)
<input type="checkbox"/> In-step	<input type="checkbox"/> In-step
<input type="checkbox"/> Back of knee	<input type="checkbox"/> Back of knee
	<input type="checkbox"/> Heel

All measurements should be in centimeters.
 * Design Pressure
 ** If measuring is done in lying position, cA please apply 0 tension
 ***cD/cG 0 tension with silicone band and straight ending
 Lateral rise: standard is 4cm AD and required is 7cm AG
 For additional product order forms, please go to:
<https://eshop.jobst-usa.com>

Comments: _____

Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application

Size	Shoe size	BNR	UOM / Box	Order Qty
X-Small	≤ 2	7965803	1	
Small	2.5 - 5.5	7965804	1	
Medium	6 - 8	7965802	1	
Large	8.5 - 11	7965902	1	
X-Large	≥ 11.5	7966001	1	



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	≤ 7.5	7957100	1	
Large	8 - 12	7957101	1	
X-Large*	≥ 13	7957102	1	



Arion Magnide® 2in1

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	< 7.5	7957110	1	
Large	≥ 8 - 10	7957111	1	
X-Large*	≥ 10.5	7957112	1	




* XL also suitable for extra wide feet with cA up to 28 cm

Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product

Size	Shoe size	BNR	UOM / Box	Order Qty
Small	< 3 - 6	7063301	1	
Medium	6.5 - 8	7063302	1	
Large	8.5 - 10.5	7063303	1	
X-Large	≥ 11	7063304	1	



Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children

Size	Age	Length	BNR	UOM / Box	Order Qty
X-Small	0-3 yrs	11.6" (29.5cm)	7957108	1	
Small	3-6 yrs	15.5" (39.5cm)	7957109	1	



Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®

Size	Circumference of lower leg	BNR	UOM / Box	Order Qty
Medium	≤ 15.7" (40cm)	7957106	1	
Large	≤ 23.6" (60cm)	7957107	1	
Extra large	≤ 39.5" (100cm)	7957107	1	



Dycem® Non-Slip Material*

- Improves grip
- Provides stability
- Helps with easier donning of a JOBST® garment

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
One size	7763600	1	



Dycem® is a registered trademark of Dycem Ltd.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
Small	131202	1 pair	
Medium	131203	1 pair	



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box	Order Qty
Standard	110913	1 donner	
Large	110910	1 donner	

